



# CAGE

## 5'-RNA sequence service/kit

New approach for genome-wide promoter identification and gene expression profiling

Quantify all RNA Polymerase II transcripts:

- ✓ Analyze transcripts not only as genes but also as Transcription Start Sites
- ✓ Discover alternative promoters for all endogenous genes
- ✓ Predict transcription factor binding motifs more reliably
- ✓ Provide a powerful tool to reveal gene network dynamics

Extensively used in research projects:

- The ENCODE project at NIH
- The FANTOM project at RIKEN
- The national project "Transcription Network Analysis" in Japan

Do you dream of a way to analyze gene expression networks?  
Discover DNAFORM and find out more about CAGE-seq!!

<https://cage-seq.com/>

# CAGE: A unique transcriptome analysis method for 5'-ends of nc/mRNA

Cap Analysis of Gene Expression (CAGE) is a new method for expression profiling and promoter identification, which allows transcriptional network analysis and transcriptome characterization.

CAGE utilizes “cap-trapping” technology to capture the 5' cap of mRNA/ncRNA. Through high volume parallel sequencing of cDNA corresponding to 5'-end of RNA and analysis of the sequenced tags, Transcription Start Sites (TSS) and transcript amount are inferred on a genome-wide scale (Fig. 1).

The process of library preparation of CAGE use neither PCR nor fragmentation of RNA which cause biased results of the gene expression in RNA-seq.




Figure 1. Overview of CAGE

## A powerful tool in gene regulation research

CAGE accurately identifies the position and expression level of TSS on a genome-wide scale with high reproducibility (Fig. 2).

The cap-trapping technology which captures one transcript with a 5' cap as one read, and the PCR-free library preparation process without fragmentation allow for digital quantification of RNA transcript abundances including low expressed transcripts.

CAGE detects precise position of TSS, which is extremely difficult for most transcriptomic technologies to identify (Fig. 3, and Table 1).

Using information obtained by CAGE, such as precise positions, variants and transcript abundances of TSS, you can detect the corresponding promoters and accurately estimate signaling cascades.




Figure 2. Correlation between replicates of CAGE  
Two CAGE technical replicates of the same sample show a very high correlation (Spearman correlation:  $R^2=0.999$ ).




Figure 3. Expression level (tags per million) of the human MYB gene of CAGE (above) and RNA-seq (below)

Table 1. Features of CAGE and other transcriptome analyses (“N/A” means not applicable)

Purpose of the study	CAGE	RNA-seq	ChIP-seq	Micro array	SAGE
<i>de novo</i> Gene/ncRNA finding	good	good	average	N/A	good
Gene expression quantification	superior	good	N/A	average	good
Determining a promoter site	superior	average	good	N/A	N/A
Motif finding for the transcription factor binding site	superior	average	superior	average	average
Identification of bidirectional enhancer RNA	superior	N/A	N/A	N/A	N/A
Determining transcription start /1 <sup>st</sup> exon site	superior	average	N/A	N/A	N/A
Determining gene structure (intron/exon, fusion gene, alternative splicing variants)	N/A	average	N/A	N/A	N/A

# CAGE identifies complex gene expression networks

CAGE accurately provides information such as alternative TSS usage on a genome-wide scale (Fig. 4). CAGE is a powerful tool in gene expression network research.

## Comparative analysis of promoter activity

Reliable estimation of promoter regions and their activities based on precise TSS information enable you to find alternative promoters and to compare promoter utilization patterns among different organs, developmental stages and diseased/normal organs which are essential in identifying expression networks underlying differentiation, development and diseases. CAGE also allows you to characterize type and genome-wide distribution of promoters (Carninci et al. 2006, Ohmiya et al. 2014).

## Valuable tool in the development of new biomarkers of cancers and other diseases

CAGE detects TSS variants of mRNAs/ncRNAs, which vary in expression level and pattern depending on the type of cancer cells, diseased/normal organs (Fig. 5). TSS variants are valuable candidate of biomarkers even in the case that there are no difference at the transcript level.

## Explore transcription factor binding motifs


CAGE enables you to explore transcription factor binding motifs. With CAGE, you can perform a genome-wide motif search around precise TSS positions which have different expression profiles depending on case and control (i.e. “up-regulated” or “down-regulated”) to obtain a list of candidates for transcription factor binding motifs which correspond to different expressions (Table 2).

The precise distance between each motif and TSS obtained by CAGE can verify estimations of associated transcription factors and enable you to find candidate genes that are regulated by specific transcription factors (Fig. 6).

## Non-coding RNA analysis


CAGE can be applied to detect and quantify long non-coding RNAs (lncRNAs), even those which are not polyadenylated. CAGE provides accurate information of the promoters and 5'-end sequence of the lncRNAs, which are difficult to determine by RNA-seq.

CAGE also enables you to identify active enhancers by detection of TSS of enhancer associated bidirectional RNAs (eRNAs). Enhancer identification by eRNAs is more precise than mapping analysis by ChIP-seq.



**Figure 4. Expression patterns of TSS variants of the P2RY6 gene between HeLa cells (cervical carcinoma cell line) and THP-1 cells (leukemia monocytic cell line)**

CAGE detects differences in the expression pattern between TSS at the single nucleotide level. Detection at this level is impossible by RNA-seq.




**Figure 5. CAGE expression pattern of a histone H3 methyltransferase gene of human CD14(+) monocytes (above) and THP-1 leukemia monocytic cells (below)**

Motif No.	Consensus Motif	Foreground: /100	Background: /1000	P value	Known Motifs (P value)
AMD_001	CAACTNGCG	27	51	1.42E-04	NA
AMD_002	GTARCNNWNSSCG	31	54	1.32E-05	NA
AMD_003	CTTCARNNNNCGA	36	108	4.77E-03	NA
AMD_004	ACGTNNNNNGNACC	28	44	1.24E-05	PPARy (9.21e-05)

**Table 2. Motif search near TSS with an expression level higher in preadipocytes than in mesenchymal stem cells**

The known adipose differentiation marker, PPARy, is detected.



**Figure 6. Kruppel-like factor 4 (KLF4)-like binding motif discovered by a motif search based on CAGE (A, below) and distance distribution of motifs from TSS (B)**

In addition to a similar sequence, the positioning of the motif relative to TSS (5B, 100bp upstream of TSS) strongly suggests that KLF4 is involved in transcription of the motif because the motif position is consistent with the reported transcription pattern of KLF4. (Figure courtesy of Ohmiya et al. *BMC Genomics*. 2014 Apr 25;15:269)

# How many TSS are present for each gene?

Riken and the FANTOM consortium have identified more than 201,802 human TSS (Science 347, 2015), which is almost four times the number of known protein-coding genes and ncRNAs, using CAGE. More than half of all known genes are regulated by multiple alternative promoters. It is strongly indicated that many of these multiple promoters affect a tissue in a specific way and may be linked to specific diseases.

CAGE peaks within 500bp of annotated 5'-end		Human		Mouse	
		Peaks	Peaks /gene	Peaks	Peaks /gene
Robust Promoter	Coding + ncRNA	184,827		116,277	
	Coding RNA	82,150	4.3	61,134	3.2
Permissive Promoter	Coding + ncRNA	1,048,124		652,860	
	Coding RNA	245,514	11.8	146,185	7.1

Andersson R, et al. *Nature* 2014 Mar 27;507(7493):455-61.

## Representative papers using CAGE

- Promoter-level transcriptome in primary lesions of endometrial cancer identified biomarkers associated with lymph node metastasis.**  
Yoshida E, et al. *Sci. Rep.* 2017 Oct 26;7:14160. doi: 10.1038/s41598-017-14418-5.
- An atlas of human long non-coding RNAs with accurate 5' ends.**  
Hon CC, et al. *Nature.* 2017 Mar 9;543(7644):199-204. doi: 10.1038/nature21374.
- Single-Nucleotide Resolution Mapping of Hepatitis B Virus Promoters in Infected Human Livers and Hepatocellular Carcinoma.**  
Altinel K, et al. *J Virol.* 2016 Nov 14;90(23):10811-10822.
- Reduced expression of APC-1B but not APC-1A by the deletion of promoter 1B is responsible for familial adenomatous polyposis.**  
Yamaguchi K, et al. *Sci. Rep.* 2016 May 24;6:26011. doi: 10.1038/srep26011.
- Enhanced Identification of Transcriptional Enhancers Provides Mechanistic Insights into Diseases.**  
Murakawa Y, et al. *Trends Genet.* 2016 Feb;32(2):76-88. doi: 10.1016/j.tig.2015.11.004.
- DeepCAGE Transcriptomics Reveal an Important Role of the Transcription Factor MAFB in the Lymphatic Endothelium.**  
Dieterich LC, et al. *Cell Rep.* 2015 Nov 17;13(7):1493-504.
- Nuclear transcriptome profiling of induced pluripotent stem cells and embryonic stem cells identify non-coding loci resistant to reprogramming.**  
Fort A, et al. *Cell Cycle.* 2015;14(8):1148-55. doi: 10.4161/15384101.2014.988031.
- Transcribed enhancers lead waves of coordinated transcription in transitioning mammalian cells.**  
Arner E, et al. *Science* 2015 Feb 27;347(6225):1010-4. doi:10.1126/science.1259418.
- Deep transcriptome profiling of mammalian stem cells supports a regulatory role for retrotransposons in pluripotency maintenance**  
Fort A, et al. *Nat. Genet.* 2014 Jun 28;46:558-566. doi:10.1038/ng.2965.
- Two independent transcription initiation codes overlap on vertebrate core promoters.**  
Haberle V, et al. *Nature* 2014 Mar 20; 507(7492):381-385. doi: 10.1038/nature12974.
- Tiny RNAs associated with transcription start sites in animals.**  
Taft RJ, et al. *Nat. Gen.* 2009 Apr 19;41(5): 572-578. doi:10.1038/ng.312.

\*For other papers, please check our website: <https://cage-seq.com/index.html>

## Project workflow

1. Total RNA submission		2. Library preparation/ Sequences		3. Bioinformatics analysis	
Amount	> 3 µg	Sequencing instrument	HiSeq 2500/ NextSeq500	Raw sequencing data (FASTQ format)	
Preferable Concentration	> 0.1 µg/µl	Amount of data	10-15 million reads/sample	Genome mapping data (BAM format)	
Purity	A260/A280 > 1.8	Sequencing	50bp/75bp Single-end	TSS cluster with annotation	
	A260/A230 > 1.8			Differential expression analysis of TSS clusters (ScatterPlot, Heat map, Clustering)	
RIN	> 7.0	<b>Total 4-6 weeks</b>		Motif search analysis	

## Ordering information

CAGE library preparation & analysis services	
Services	Price <sup>1</sup>
Library preparation <sup>2</sup>	500 USD/sample
Sequencing	250 USD/sample
Bioinformatics analysis	250 USD/sample

1. Shipping : 200 USD/shipment

2. Library is prepared for Illumina-platform

CAGE library preparation Kit <sup>3</sup>		
Package size	Price	Cat. No.
8 samples	2,000 USD	52003-8
48 samples	10,000 USD	52003-48

3. Shipping : 800 USD/shipment

4. Library is prepared for Illumina-platform



**K.K. DNAFORM**

Life science research center 402, 1-1-43 Suehiro-cho, Tsurumi-ku, Yokohama, Kanagawa 230-0045 Japan  
TEL: +81-45-508-1539 FAX: +81-45-510-0608 E-mail: [contact@dnaform.jp](mailto:contact@dnaform.jp)

