

Quantify all mRNA / ncRNA during transcription New methods for gene expression

Analysis to enable identification of active enhancers

- ✓ Novel gene expression analysis based on nascent RNA (NET-RNA)
- Captures short-lived eRNA and identifies enhancer regions
- Gene expression analysis based on instantaneous transcriptional activity
- Just send the cells or tissues to us, delivered with Bioinformatics analysis



Discovery of short-lived non-coding RNAs

In addition to mRNA, RNA polymerase II is known to transcribe other RNAs (long-non-coding RNAs) such as enhancer RNAs (eRNAs), which are expressed in both directions from the enhancer region, and promoter upstream transcripts (PROMPTs), which are transcribed in the opposite direction to mRNA from upstream of the promoter. eRNAs and PROMPTs have a 5'-end cap structure, are extremely short-lived and are degraded in the nucleus immediately after synthesis.



NET-CAGE : New transcriptome analysis to capture nascent strand RNA

NET-CAGE (Native Elongating Transcript Cap Analysis of Gene Expression)¹ is a technology that concentrates nascent RNA in the nucleus and maps the position of transcription start sites of the nascent RNA on the genome with accuracy of single nucleotide resolution. Not only mRNA, but also short-lived RNAs such as eRNAs and PROMPTs can be captured and quantified.



1: Hirabayashi et al. Nature Genetics 2019 volume 51, 1369-1379

What is CAGE?

CAGE (Cap Analysis of Gene Expression) is a technology for determine sequence of the 5' end of RNA and mapping it onto the genome based on the "Cap-trapping" method, which captures the Cap structure of RNA polymerase II transcripts.



Identify the transcription start sites and quantify each transcript expression levels!

Detection of active enhancer and promoter regions by NET-CAGE

Transcription start sites of eRNAs that are expressed in both directions from enhancer regions, enables the identification of active enhancers in high resolution. PROMPTs can also be used to identify the extent of activated promoter regions. Furthermore, NET-CAGE capture mRNAs in transcription at the same time, enables to quantify based on the true transcript expression levels of mRNAs.



Captures eRNAs and PROMPTs, detects activated enhancer and promoter regions on a genome-wide scale!





High-resolution identification of individual enhancer regions in known super-enhancer by NET-CAGE!



Two replicates of NET-CAGE were performed on the same sample and the expression level were plotted(Log_2CPM).High correlation coefficients are shown.

NET-CAGE Features

- 1 Activated enhancer and promoter regions can be identified at high resolution.
- 2 Unlike GRO-seq and PRO-seq, analysis is possible from frozen samples and tissue samples. * 1
- 3 Gene expression comparison data based on real-time transcriptional activity is also obtain at the same time.
- 4 Just send the cells or tissue to us. Delivered with Bioinformatics analysis.

* 1 Some tissue types, such as bone and cartilage, may be difficult to analyze with NET-CAGE. Please consult in advance.

Project workflow

Table 1. Features of NET-CAGE and other transcriptome analyses ("N/A" means not applicable)						
Purpose of the study	NET- CAGE	CAGE	RNA- seq	ChIP- seq		
de novo Gene / ncRNA finding	good	good	good	average		
Gene expression quantification	superior	superior	good	N/A		
Determing a promoter sites	superior	superior	average	good		
Motif finding for the transcription factor binding site	superior	superior	average	superior		
Determing the eRNA / PROMPTs	superior	average	N/A	N/A		
Determing transcription start / 1 st exon site	superior	superior	average	N/A		
Determing gene structure (intron/exon, fusion gene, alternative splicing variants)	N/A	N/A	average	N/A		
Difficulty of sample adjustment	average	average	good	N/A		
Gene expression comparison based on instantaneous transcriptional activity	superior	N/A	N/A	N/A		

1. Sample submision		2. Library preparation/ sequences		3. Bioinfomatics analysis	
Sample Cells or tissues		NGS	NextSeq2500/NextSeq500	Raw sequencing data (FASTQ format)	
Number of cells (for cultured cells)	About 2×107 cell*2	Amount of data	10-15million reads/sample	Genome mapping data (BAM format)	
Sample condition	Frozen at -80°C	Sequencing	50bp/75bp Single-end	TSS cluster with anotation	
*1 Some tissues may not be ar Please consult us in advanc	nalyzed. The together with the sample	e volume.		Differential expression analysis of TSS Clusters ScatterPlot, Heat map, Clustering)	
*2 If possible, please send two	tubes of 2×10^7 cells.	Total	5-8 weeks	Comparison of expression with enhancers detected in the RIKEN FANTOM project. Motif search analysis	

Ordering information

CAGE library preparation & analysis services				
Services	Price ¹			
NET-RNA extraction	150 USD/sample			
Library preparation ²	500 USD/sample			
Sequencing	250 USD/sample			
Bioinformatics analysis	250 USD/sample			

1. Shipping : 200 USD/shipment

2. Library is prepared for Illumina-platform

YouTube Check out our YouTube channel! **Technical webinar about NET-CAGE**



https://www.youtube.com/@user-ky4xw1tz2i/videos





K.K. DNAFORM

Ask Sanshin Building 3F, 2-6-29 Tsurumi-chuo, Tsurumi-ku, Yokohama, Kanagawa 230-0051 Japan TEL: +81-45-508-1539 FAX: +81-45-510-0608 E-mail: contact@dnaform.jp

